

somnoware

Physician Quickstart Guide

Why Project USER?

- Move platform off of Flash technology to prepare for removal of support by browsers on 1/1/2021
- Improve the platform by making it more User and task centric, more Secure and drive more Efficiency


Consistent Patterns

In addition to starting with user tasks to bring work to one place, we made every effort to use consistent patterns to make the platform easy to learn and use. These patterns will be noted as we way through the core workflow for reading and signing studies.

1 - Login – App3.Somnoware.com

Use your current credentials (user name and password.) In near term, you can log into today's platform (App.Somnoware.com) or future App3 without issue. They work off the same data Note: you will be taken to Physician Study Manager upon login

2 - Physician Study Manager Orientation



2.1 - Left side main navigation in dark blue vs. top - "Patient", "Schedule", "Study Manager", Etc.


2.2 - Sub-navigation on left in white is like old filters. Below seeing patients Pending Interpretation

2.3 - Facility, logout, search, "My Studies" and Date range - Upper right



3 - Primary and secondary actions in study manager

They are viewable when hovering over a patient and secondary are accessible by clicking the three dots on far right.



3.1 - Primary actions


Patient Profile
Study notes
Interpretation (= View report)
View tracking (= Cloud PSG).

3.2 - Secondary actions

Depends on the study status.

Note: no double clicks, only single clicks in the new platform

4 - Reading & Signing a study – click on “Interpretation” per (3) above.



4.1 - Upper right can put in “Edit” mode.

Editor is the same as in the current Flex app.

4.2 - Lower right to “Sign.”

4.3 - Across top can tab to “Scored report”, “Tech note”, “Documents.”

4.4 - Upper right you can open tracing in new window with “View tracings.”

If you have any feedback or need some help please email use at

Betahelp@somnoware.com

Thanks for your time. The Somnoware Team